

Safe automatic disinfection

Manual disinfection is tedious and always leaves room for error.
Automate your disinfection to **save time and ensure total safety**.
The Mobiwatch device and product combination disinfects all your accessible and inaccessible surfaces, as well as all equipment within the treated area.

Maxibio solution complies with European standard EN 17 272

The Mobiwatch Maxibio diffusion system

Diffuses a 'dry mist' into the atmosphere that saturates the space and eradicates pathogenic microorganisms.

Disinfects all surfaces, both accessible and inaccessible!



2-year warranty

Features

- Compact unit
- Fully automatic
- No need to tie up staff
- Suitable for all types of premises
- Optimised amount of disinfectant dispensed
- Rapid fogging
- Disinfection traceability

Technical Specifications	Maxibio
Dispersion system	Venturi effect – suction
Tank capacity	2 L
Volume to be treated	5 to 320 m ³ with Mobiwatch® SPP 5 to 165 m ³ with O2SAFE 7.4®
Minimum diffusion rate	1800 mL/hour
Compressor – Nozzle (316 stainless steel)	Robust and simple system
Power supply	230 V – 50/60 Hz – 310 W
Dimensions (mm) – Weight	345 x 325 x 345 – 12 kg

Mobiwatch® SPP is a disinfectant based on hydrogen peroxide and peracetic acid.

The ready-to-use Mobiwatch SPP disinfectant features a unique, patented formulation. 5L containers.

Mobiwatch SPP: < 8% H₂O₂ + < 0.2% APA

Mobiwatch SPP is effective against all four spectra
BACTERICIDAL – VIRUCIDAL – FUNGICIDAL – SPORICIDAL



2-year shelf life

